

Protein serine/threonine kinase

Chemical Properties

CAS No. : 9026-43-1

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Protein serine/threonine kinase (Ser/Thr protein kinase) is a biochemical reagent that serves as a biomaterial or organic compound, utilized in life science research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481