

## Fmoc-Cys-Asp10 TFA

## Chemical Properties

CAS No. :

Formula: C60H68F3N11O36S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Fmoc-Cys-Asp10 (TFA) is a non-releasable oligopeptide linker involved in synthesizing releasable oligopeptide linkers. These releasable linkers are used to deliver drugs to fracture-targeting oligopeptides, thereby reducing the healing time of fractured femurs.
Targets(IC50)	Amino Acids and Derivatives

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481