

[Ala18] Endothelin-1, human

Chemical Properties

CAS No. :

Formula: C108H163N25O30S5

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	[Ala18] Endothelin-1, human, is a bioactive peptide utilized in various biochemical research applications.
Targets(IC50)	Endothelin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481