


Agouti-Related Protein (25-51) (human)

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 209460-61-7 |
| Formula: | C130H221N37O35S |
| Molecular Weight: | 2894.44 |
| Storage: | Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small> |



Biological Description

| | |
|---------------|--|
| Description | Agouti-Related Protein (25-51) (human) (AGRP (25-51) (human)) corresponds to the amino-terminal fragment comprising amino acids 25 to 51 of the agouti-related protein (AGRP). |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 0.3455 mL | 1.7274 mL | 3.4549 mL |
| 5 mM | 0.0691 mL | 0.3455 mL | 0.691 mL |
| 10 mM | 0.0345 mL | 0.1727 mL | 0.3455 mL |
| 50 mM | 0.0069 mL | 0.0345 mL | 0.0691 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481