

## Pepinh-Control

## Chemical Properties

CAS No. :

Formula: C172H274N54O38S2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Pepinh-Control is a negative control peptide.
Targets(IC50)	Others
In vitro	Pepinh-Control (100 $\mu$ M, 12 h) serves as a negative control for Pepinh-MYD in the pre-cultivation of HDPFs (human dental pulp fibroblasts).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481