

## Minigastrin

## Chemical Properties

CAS No. :	60748-07-4
Formula:	C74H99N15O26S
Molecular Weight:	
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Minigastrin is an acidic peptide that is useful for studying peptide-peptide interactions.
Targets(IC50)	Cholecystokinin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481