

VPSTPPTPSPSTPPTSPS TFA

Chemical Properties

CAS No. :

Formula: C₈₁H₁₂₈N₂₀O₂₈.xC₂H₅F₃O₂

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	VPSTPPTPSPSTPPTSPS (TFA) is a peptide from the IgA1 hinge region.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481