

FRET5-VWF73

Chemical Properties

CAS No. :

Formula: C370H583N103O113S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FRET5-VWF73 is a peptide composed of 73 amino acids, serving as a fluorescent substrate for detecting ADAMTS13 [Ex=340 nm; Em=450 nm].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481