

5-FAM-PMDM6

Chemical Properties

CAS No. :

Formula: C84H96ClN14O24PS

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	5-FAM-PMDM6 (PMDM6-F) is a biologically active peptide. It serves as a fluorescently labeled probe used in MDM2 binding assays.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481