

A β 1-14- ϵ K-HBsAg3 Th

Chemical Properties

CAS No. :

Formula: C175H287N49O51

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	A β 1-14- ϵ K-HBsAg3 Th (UB311 immunogen II) is a bioactive peptide. It is one of the A β 1-14 targeting peptides (B cell epitope) for UB-311.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481