

AChR $\alpha$ (97-116) TFA

## Chemical Properties

CAS No. :

Formula: C106H162N24O30.xC2HF3O2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AChR $\alpha$ (97-116) TFA is a polypeptide useful in inducing experimental autoimmune myasthenia gravis (EAMG).
Targets(IC50)	AChR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481