

$\alpha\beta 6$ -BP TFA

Chemical Properties

CAS No. :

Formula: C₉₃H₁₆₄N₃₄O₂₇.xC₂H₅F₃O₂

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	$\alpha\beta 6$ -BP TFA is a selective $\alpha\beta 6$ -binding peptide used for molecular imaging research.
Targets(IC ₅₀)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481