

m-PEGn-NHS ester

Chemical Properties

CAS No. :	174569-25-6
Formula:	(C ₂ H ₄ O) _n C ₈ H ₁₁ N ₅
Molecular Weight:	1787.086
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	m-PEGn-NHS ester is a PROTAC linker, belonging to the PEG category, and is utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5596 mL	2.7978 mL	5.5957 mL
5 mM	0.1119 mL	0.5596 mL	1.1191 mL
10 mM	0.056 mL	0.2798 mL	0.5596 mL
50 mM	0.0112 mL	0.056 mL	0.1119 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481