

## Diethylene glycol diacetate

## Chemical Properties

CAS No. : 628-68-2


Formula: C<sub>8</sub>H<sub>14</sub>O<sub>5</sub>

Molecular Weight: 190.19

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Diethylene glycol diacetate is a PROTAC linker classified under the PEG category, utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2579 mL	26.2895 mL	52.579 mL
5 mM	1.0516 mL	5.2579 mL	10.5158 mL
10 mM	0.5258 mL	2.6289 mL	5.2579 mL
50 mM	0.1052 mL	0.5258 mL	1.0516 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481