

TCO4-PEG3-Maleimide

Chemical Properties

CAS No. : 1914971-04-2


Formula: C₂₄H₃₇N₃O₈

Molecular Weight: 495.57

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TCO4-PEG3-Maleimide is a PROTAC linker of the PEG class. It contains TCO and Maleimide groups, which can undergo specific "click" reactions with tetrazine groups or thiol groups, respectively, or participate in "thiol-Michael" addition.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0179 mL	10.0894 mL	20.1788 mL
5 mM	0.4036 mL	2.0179 mL	4.0358 mL
10 mM	0.2018 mL	1.0089 mL	2.0179 mL
50 mM	0.0404 mL	0.2018 mL	0.4036 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481