

## (R)-TCO4-PEG3-Maleimide

### Chemical Properties

CAS No. :

Formula: C<sub>24</sub>H<sub>37</sub>N<sub>3</sub>O<sub>8</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	TCO4-PEG3-Maleimide is a PROTAC linker from the PEG category. It incorporates TCO and Maleimide groups, which can selectively undergo "click" reactions with tetrazine groups or react with thiol groups through a "thiol-maleimide" reaction.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481