

CG428-NEG

Chemical Properties


CAS No. : 2714560-39-9

Formula: C44H45FN10O6

Molecular Weight: 828.89

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CG428-NEG is a negative control for TRK degraders. CG428 inhibits cell growth and reduces levels of TPM3-TRKA.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2064 mL	6.0322 mL	12.0643 mL
5 mM	0.2413 mL	1.2064 mL	2.4129 mL
10 mM	0.1206 mL	0.6032 mL	1.2064 mL
50 mM	0.0241 mL	0.1206 mL	0.2413 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481