

Pomalidomide 4'-PEG6-azide

Chemical Properties

CAS No. : 2624181-62-8


Formula: C₂₇H₃₈N₆O₁₀

Molecular Weight: 606.62

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomalidomide 4'-PEG6-azide (compound 33e) is an E3 ligase ligand-linker conjugate designed to recruit the CRBN protein, and it can be utilized in the synthesis of PROTAC.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6485 mL	8.2424 mL	16.4848 mL
5 mM	0.3297 mL	1.6485 mL	3.297 mL
10 mM	0.1648 mL	0.8242 mL	1.6485 mL
50 mM	0.033 mL	0.1648 mL	0.3297 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481