


Boc-Leu-Gly-Arg-pNA acetate

Chemical Properties

CAS No. :	176786-91-7
Formula:	C ₂₇ H ₄₄ N ₈ O ₉
Molecular Weight:	624.69
Storage:	Keep away from direct sunlight Store at -20°C <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Boc-Leu-Gly-Arg-pNA is a chromogenic synthetic substrate.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6008 mL	8.004 mL	16.0079 mL
5 mM	0.3202 mL	1.6008 mL	3.2016 mL
10 mM	0.1601 mL	0.8004 mL	1.6008 mL
50 mM	0.032 mL	0.1601 mL	0.3202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481