

## Sudan red 7B

## Chemical Properties

CAS No. : 6368-72-5


Formula: C<sub>24</sub>H<sub>21</sub>N<sub>5</sub>

Molecular Weight: 379.46

Storage: Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sudan Red 7B is a red non-fluorescent dye used for staining lipid bodies.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6353 mL	13.1766 mL	26.3532 mL
5 mM	0.5271 mL	2.6353 mL	5.2706 mL
10 mM	0.2635 mL	1.3177 mL	2.6353 mL
50 mM	0.0527 mL	0.2635 mL	0.5271 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481