

Distearin

Chemical Properties

CAS No. : 1323-83-7

Formula: C₃₉H₇₆O₅

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Distearin is a diacylglycerol, where stearic acid is present at two positions. Stearic acid is a long-chain dietary saturated fatty acid found in various animal and plant fats. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481