


Spinosyn A aglycone

Chemical Properties

CAS No. :	149560-97-4
Formula:	C ₂₄ H ₃₄ O ₅
Molecular Weight:	402.52
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Spinosyn A aglycone ((-)-Spinosyn A aglycon) is the glycosidic form of the insecticide Spinosyn A. It is an acidic degradation product resulting from the hydrolysis of two sugar moieties on Spinosyn A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4843 mL	12.4217 mL	24.8435 mL
5 mM	0.4969 mL	2.4843 mL	4.9687 mL
10 mM	0.2484 mL	1.2422 mL	2.4843 mL
50 mM	0.0497 mL	0.2484 mL	0.4969 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481