


## Salvinorin B

## Chemical Properties

CAS No. : 92545-30-7  
 Formula: C<sub>21</sub>H<sub>26</sub>O<sub>7</sub>  
 Molecular Weight: 390.43  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Salvinorin B is a metabolite produced by <i>Salvia divinorum</i> and is also the primary metabolite of Salvinorin A in the blood of rhesus monkeys.
Targets(IC50)	Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5613 mL	12.8064 mL	25.6128 mL
5 mM	0.5123 mL	2.5613 mL	5.1226 mL
10 mM	0.2561 mL	1.2806 mL	2.5613 mL
50 mM	0.0512 mL	0.2561 mL	0.5123 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481