

ALK/ROS1-IN-3

## Chemical Properties

CAS No. :

Formula: C32H32N4O2

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | ALK/ROS1-IN-3 (compound C01) is a dual inhibitor targeting ROS1 and ALK, with IC50 values of 42.3 nM for ROS1G2032R and 49.1 nM for ALKG1202R. |
| Targets(IC50) | ALK,ROS Kinase   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481