

DA5-HTL

## Chemical Properties

CAS No. :

Formula: C<sub>39</sub>H<sub>58</sub>Cl<sub>2</sub>N<sub>8</sub>O<sub>8</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DA5-HTL is a compound that combines dasatinib with the HaloTag system, effectively inhibiting the growth of tumor cells, with a GI <sub>50</sub> of 1.923 nM.
Targets(IC <sub>50</sub> )	Apoptosis,Bcr-Abl,Autophagy,Src

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481