

ssVACV-70mer sodium

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ssVACV-70mer sodium is a 70 bp single-stranded oligonucleotide containing a viral DNA motif derived from vaccinia virus DNA. Unlike its double-stranded counterpart, dsVACV 70mer, ssVACV-70mer sodium does not act as an IFN inducer.
Targets(IC50)	STING

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481