

dsVACV-70mer sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	dsVACV-70mer sodium is a 70 base pair double-stranded oligonucleotide that contains a viral DNA motif derived from vaccinia virus DNA. It effectively induces IFN- β through a STING-dependent mechanism.
Targets(IC50)	IFNAR,STING

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481