

ORN 06 sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ORN 07 (sodium) is a synthetic single-stranded RNA rich in uracil (U) and consists of six repeating UUAU sequence motifs. AU-rich oligoribonucleotides (ORN) are capable of stimulating TLR7.
Targets(IC50)	TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481