

PDH E1-IN-1

Chemical Properties

CAS No. :

Formula: C₂₃H₂₈N₄O₂S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PDH E1-IN-1 is a selective TPP-competitive inhibitor of PDH E1 with an IC ₅₀ value of 0.99 μ M.
Targets(IC ₅₀)	Dehydrogenase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481