

$\beta$ -catenin/BCL9 PPI-IN-1

## Chemical Properties

CAS No. :

Formula: C<sub>22</sub>H<sub>16</sub>O<sub>7</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	$\beta$ -catenin/BCL9 PPI-IN-1 (compound B4) is an effective inhibitor of the interaction between $\beta$ -catenin and B-cell lymphoma 9 proteins ( $\beta$ -catenin/BCL9 PPI), with an IC <sub>50</sub> value of 2.25 $\mu$ M.
Targets(IC <sub>50</sub> )	Wnt/beta-catenin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481