

SB1617

Chemical Properties

CAS No. :

Formula: C₂₂H₂₁ClN₆O₅S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SB1617 is a neuroinflammation modulator with neuroprotective effects, reducing pathogenic tau levels through anti-inflammatory activity mediated by microglia.
Targets(IC50)	Microtubule Associated

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481