


QP5038

Chemical Properties

CAS No. : 3031429-88-3
 Formula: C₂₁H₂₀FN₅
 Molecular Weight: 361.42
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	QP5038 is an inhibitor of QPCTL, exhibiting an IC ₅₀ value of 3.8 nM. It also possesses antitumor properties.
Targets(IC ₅₀)	Beta Amyloid

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7669 mL	13.8343 mL	27.6686 mL
5 mM	0.5534 mL	2.7669 mL	5.5337 mL
10 mM	0.2767 mL	1.3834 mL	2.7669 mL
50 mM	0.0553 mL	0.2767 mL	0.5534 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481