


## FXR antagonist 3

## Chemical Properties

CAS No. : 3024652-96-5  
 Formula: C<sub>32</sub>H<sub>36</sub>N<sub>2</sub>O<sub>4</sub>  
 Molecular Weight: 512.64  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	FXR antagonist 3 (V02-8) is an intestine-specific FXR receptor antagonist with an IC <sub>50</sub> value of 0.89 μM.
Targets(IC <sub>50</sub> )	FXR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9507 mL	9.7534 mL	19.5069 mL
5 mM	0.3901 mL	1.9507 mL	3.9014 mL
10 mM	0.1951 mL	0.9753 mL	1.9507 mL
50 mM	0.039 mL	0.1951 mL	0.3901 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481