

ALF-PD-FAPI

Chemical Properties

CAS No. :

Formula: C43H54AlF3N10O16

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ALF-PD-FAPI exhibits an affinity for FAP with an IC50 value of 0.13 nM. It is characterized by specific uptake, high internalization rate, and low extracellular excretion in vitro. ALF-PD-FAPI can serve as a FAP-targeting tracer.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481