

DAT-IN-1

Chemical Properties

CAS No. :

Formula: C₂₉H₃₄F₂N₂O₂S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | DAT-IN-1 (Compound 12b) is an atypical dopamine transporter (DAT) inhibitor with a K _i of 8.37 nM. It can be utilized in research focused on psychostimulant use disorder (PSUD). |
| Targets(IC50) | Dopamine Receptor |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481