

DAT-IN-1

## Chemical Properties

CAS No. :

Formula: C<sub>29</sub>H<sub>34</sub>F<sub>2</sub>N<sub>2</sub>O<sub>2</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DAT-IN-1 (Compound 12b) is an atypical dopamine transporter (DAT) inhibitor with a K <sub>i</sub> of 8.37 nM. It can be utilized in research focused on psychostimulant use disorder (PSUD).
Targets(IC50)	Dopamine Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481