

SiR-PEG4-BCN

Chemical Properties

CAS No. :

Formula: C48H62N4O9Si

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SiR-PEG4-BCN is a fluorescent dye and a reagent for click chemistry. It includes a BCN group that facilitates strain-promoted alkyne-azide cycloaddition (SPAAC) with molecules containing azide groups.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481