

T-1-MBHEPA

## Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>25</sub>N<sub>7</sub>O<sub>4</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	T-1-MBHEPA is a potent VEGFR-2 inhibitor with an IC <sub>50</sub> of 0.121±0.051 μM. It exhibits both antitumor and anti-angiogenic properties.
Targets(IC <sub>50</sub> )	VEGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481