

RF16

## Chemical Properties

CAS No. :

Formula: C46H64BClO15

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RF16 is a fluorescent polymer material. It is precisely retained within cells and used to track the IRE1-Halo protein via the HaloTag ligand.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481