

## TAMRA-PEG4-Me-Tet

## Chemical Properties

CAS No. :

Formula: C46H52N8O9

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TAMRA-PEG4-Me-Tet is a TAMRA dye derivative that includes four PEG units. TAMRA-PEG8-Me-Tet contains a Tetrazine group, which can participate in an inverse electron demand Diels-Alder reaction (IEDDA) with molecules containing a TCO group.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481