

AF488 streptavidin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AF488 streptavidin is a fluorescently labeled streptavidin. It consists of a biotin-binding protein (streptavidin) covalently linked to a fluorescent label (AF488). AF488 is a bright and photostable green fluorophore with emission and excitation wavelengths (λ_{em} and λ_{ex}) of 520 nm and 470 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481