

AF488 carboxylic acid

Chemical Properties

CAS No. :

Formula: C₂₁H₁₁K₃N₂O₁₁S₂

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AF488 carboxylic acid represents the non-reactive form of AF 488. The emission (λ_{em}) and excitation (λ_{ex}) wavelengths for AF488 are 520 nm and 470 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481