

## AF488 carboxylic acid

### Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>11</sub>K<sub>3</sub>N<sub>2</sub>O<sub>11</sub>S<sub>2</sub>

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	AF488 carboxylic acid represents the non-reactive form of AF 488. The emission ( $\lambda_{em}$ ) and excitation ( $\lambda_{ex}$ ) wavelengths for AF488 are 520 nm and 470 nm, respectively.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481