

AF488 amine

Chemical Properties

CAS No. :

Formula: C₃₃H₄₃N₅O₁₀S₂

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AF488 amine is a bright green fluorescent dye used for labeling various molecules, such as antibodies. This compound has an emission wavelength (λ_{em}) of 520 nm and an excitation wavelength (λ_{ex}) of 470 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481