

## Preactivated APC-Cy5.5 Maleimide

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Preactivated APC-Cy5.5 Maleimide is a thiol-reactive dye that reacts with free thiol groups on proteins. Once activated, it readily conjugates to proteins or antibodies without altering the spectral characteristics of APC-Cy/YF.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481