

Preactivated PE-Cy5 Maleimide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Preactivated PE-Cy5 Maleimide is a thiol-reactive dye that can react with free thiol groups on proteins. Once activated, it readily binds to proteins or antibodies without altering the spectral characteristics of APC-Cy/YF.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481