

## Fluorescent ACKR3 antagonist 1

### Chemical Properties

CAS No. :

Formula: C48H62BF3N6O6

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Fluorescent ACKR3 antagonist 1 (compound 18a) functions as an atypical chemokine receptor 3 (chemokine receptor3) antagonist.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481