

SiR-PEG4-NHS ester

Chemical Properties

CAS No. :

Formula: C₄₂H₅₂N₄O₁₁Si

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SiR-PEG4-NHS ester is a near-infrared fluorescent dye extensively used for labeling in cell dyes, biological dyes, biomolecules, and particles. [Ex/Em = 652/674 nm]
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481