

Sulfo-Cy7 tetrazine

Chemical Properties

CAS No. :

Formula: C47H52KN7O7S2

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cy7 tetrazine is a near-infrared, water-soluble fluorescent moiety. It contains a methyl tetrazine group, allowing it to swiftly and efficiently couple with cycloalkenes through a metal-free reaction known as TCO ligation.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481