

AgAuSe-PEG-NH2 (1000 nm)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	AgAuSe-PEG-NH2 (1000 nm) is a fluorescent quantum dot that emits in the near-infrared II region, with an emission peak at 1000 nm, suitable for bioimaging applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481